

Buffer, Inoculum and Sample Preparation

07/15/2014

RF

Page 1 of 3

Buffer Preparation: Kansas State Buffer (see ref. 5)

Solution A	grams / liter
KH ₂ PO ₄	10.0
MgSO ₄ ·7 H ₂ O	0.5
NaCl	0.5
CaCl ₂ ·2 H ₂ O	0.1
Urea (optional)	0.5

Solution B	grams / 100 ml
Na ₂ CO ₃	15.0
Na ₂ S·9 H ₂ O	1.0

Mix 20 ml of Solution B with 1000 ml of Solution A and adjust pH to 6.8 by adding Solution B before each use.

Buffer Preparation: Goering - Van Soest Buffer (see ref. 3)**In vitro rumen buffer solution**

Distilled water	1 liter
NH ₄ HCO ₃	4 g
NaHCO ₃	35 g
Resaruzin 0.1% (w/v)	
Dissolve 0.1 g resaruzin into 100 ml dH ₂ O	

In vitro rumen macromineral solution

Distilled water	1 liter
Na ₂ HPO ₄ anhydrous	5.7 g
KH ₂ PO ₄ anhydrous	6.2 g
MgSO ₄ ·7 H ₂ O	.6 g

In vitro micromineral solution

CaCl ₂ ·2 H ₂ O	13.2 g
MnCl ₂ ·4H ₂ O	10.0 g
CoCl ₂ ·6 H ₂ O	1.0 g
FeCl ₃ ·6 H ₂ O	8.0 g

Reducing solution

Cysteine solution	625.0 mg
1N NaOH	4.0 ml
Distilled Water	95 ml
Na ₂ S·9H ₂ O	625.0 mg

Add to volumetric and bring volume to 100 ml with distilled water

Buffer Preparation: Cone's Buffer / Mineral Solution (see ref. 4)

	grams / liter
NaHCO ₃	8.75
NH ₄ HCO ₃	1.00
Na ₂ HPO ₄	1.43
KH ₂ PO ₄	1.55
MgSO ₄ ·7H ₂ O	0.15
Na ₂ S	0.52
CaCl ₂ ·2 H ₂ O	0.017
MnCl ₂ ·4 H ₂ O	0.015
CoCl ₂ ·6 H ₂ O	0.002
FeCl ₃ ·6 H ₂ O	0.012
Resazurin (optional)	0.125

Buffer, Inoculum, and Sample Preparation

For Systems purchased after 7/09 (Modules will have RED caps)

Buffer Preparation Procedure

- (1) Maintain all glassware at 39°C.
- (2) Mix enough buffer solution to support planned number of Gas Production Modules. Each Module's Glass Bottle can contain up to 200 ml of solution. (If using the larger 500 ml or 1000 ml Glass Bottle, make calculations using appropriate ratios). A smaller amount may be used to provide larger head space and accommodate a larger sample size.
- (3) Prepare and pre-warm your buffer solutions. To analyze 20 samples using the buffer formula provided, combine ~340 ml of solution B to 1600 ml of solution A (1:5 ratio). The exact amount of A to B should be adjusted to obtain a final pH of 6.8 at 39°C. No further adjustment of pH is necessary.
- (4) Add 100 ml of combined pH corrected solution to each Module's Glass Bottle and allow temperature of Glass Bottle and buffer to equilibrate for 20 to 30 minutes at 39°C.
- (5) Add sample to Module's Glass Bottle and allow to equilibrate with buffer solution.
- (6) While the buffer is equilibrating prepare the rumen inoculums.

Inoculum Preparation Procedure

- (1) Preheat two 2L thermos bottles by filling with 39°C water.
- (2) Empty heated water just prior to collection of rumen inoculum.
- (3) Using the appropriate collection procedure, remove 600 to 1,000 ml of rumen inoculum and place in thermos. Include approximately two "fistfuls" of the fibrous mat from the rumen with your collection in one thermos.
- (4) Empty the rumen inoculum and fibrous mat from the thermoses into a blender.
- (5) Purge the blender container with CO₂ gas and blend at a high speed for 30 seconds. The blending action serves to dislodge microbes that are attached to the mat and assure a representative microbial population for the fermentation.
- (6) Filter the blended digesta through 4 layers of cheesecloth into a pre-heated (39°C) flask. NOTE: Allow for extra cheesecloth around the edges to facilitate squeezing contents of filtered mat.
- (7) The flask should be continually purged with CO₂ before and after the transfer of the inoculum.
- (8) Be careful to maintain temperature.
- (9) Add the 25 ml of inoculum to the equilibrated buffer solution and sample in each Module's Glass Bottle.
- (10) Purge the Glass Bottle with CO₂ gas for 30 seconds. A purge system can be purchased from ANKOM Technology which allows the purging to be accomplished in a closed Module. If you don't have an ANKOM purge system then remove the Glass Bottle and allow CO₂ to flow into it.
- (11) Repeat process for all Glass Bottles to be used.

NOTE:

Do not allow CO₂ gas to bubble through the buffered inoculum. Instead, use the CO₂ to form a gaseous blanket over the contents of the Glass Bottle.

Buffer, Inoculum, and Sample Preparation

For Systems purchased after 7/09 (Modules will have RED caps)

Sample Preparation

The quantity of sample (substrate) to add to the Module will vary from 1g for a forage sample with minimal fermentable substrate to .25 g for a highly fermentable substrate such as a high starch diet. In addition to the quantity of fermentable substrate, the length of time the incubation is conducted must be taken into consideration. The quantities of substrate and buffer may be sufficient for 24 hour incubation but may not be sufficient for a 48 hour incubation. By measuring the pH at the end of the incubation period you can determine if the buffer maintained the proper pH throughout the incubation. This will allow you to alter the substrate-to-sample ration to fit the desired incubation period.

With the samples, buffer and rumen inoculum in place, be sure the digestion jars (Glass Bottles) are maintained at temperature (39o C) over the time desired for the data capture.

References

1. Theodorou, M.K., Lowman, R.S., Davies, Z.S., Cuddelford, D., and Owen. E., 1998, Principles of techniques that rely on gas measurement in ruminant nutrition. Occasional Publication No. 22, British Society of Animal Science. p. 55.
2. Pell, A.N., Pitt, R.E., Doane, P.H., and Schofield, P., 1998, The development, use and application of gas production technique at Cornell University, USA, p.45.
3. Goering, H.K. and VanSoest, P.J., 1970, Forage fiber analysis (apparatus, reagents, procedures and some applications), Agricultural Handbook No. 379 ARS-USDA, Washington, DC.
4. Cone, J.W., 1998, The development, use and application of the gas production technique at the DLO Institute for Animal Science and Health (IO-DOL), Lelystad, The Netherlands, Occasional Publication No. 22, British Society of Animal Science. p. 65.
5. Marten, G.C. and Barnes, R.F., 1980, Prediction of Energy Digestibility of Forages with In Vitro Rumen Fermentation and Fungal Enzyme Systems, in Standardization of analytical methodology for feeds: Proceedings of a workshop held in Ottawa, Canada. 12-14 March 1979. Ottawa, Ont. IDRC.