




1. Prepare the instrument for service.
 - A. Power off the instrument and unplug the instrument.
 - B. Remove the back and top panels of the electrical enclosure by removing the 14 screws.
 - C. Once open you will find the MUX-T board on the side of the electrical enclosure marked with the blue arrow.
 - D. Use a grounded anti-static wrist strap to prevent static damage to the circuit board.

2. Identify the MUX-T board. Disconnect the J2, J3, J4, J5, J7, J8 and J10 (3-pin white connectors). Remove the 4-pin connector near the top of the board. Remove the power leads labeled J-1T. Also remove the 2-pin connector near the bottom of the board.



3. Unscrew the 4 nuts securing the board to the side panel of the electrical enclosure. Be sure to save the 4 nuts and 4 plastic stand-offs. Replace with the new MUX-T board. Secure board with the stand-offs behind the board and 4 nuts.
4. Reattach the wiring connectors to the board as shown.



5. From the Touch Screen Display, go to Diagnostics / Temperatures. Check the temperatures to be sure the individual paddle temperature readings and the in-line heater are reporting temperature correctly. The screw on the blue box (upper right corner of MUX-T) is the adjustment screw for the in-line heater temperature. Adjust this to match the paddle temperatures.
6. Close up the electrical enclosure with top and back panels using the 14 screws removed earlier.

