

This Service Procedure assists the user and ANKOM Technical Service diagnose issues that may arise with the ANKOM XT10 and XT15 Fat Analyzers.

Customer Name: \_\_\_\_\_ Date: \_\_\_\_\_

Location: \_\_\_\_\_ E-Mail: \_\_\_\_\_

Instrument Serial #: \_\_\_\_\_ Best phone # to reach you: \_\_\_\_\_

1. What version of the program is on your instrument? \_\_\_\_\_  
(It is displayed on the LCD screen immediately after turning the instrument on, after the instrument name, XT10 or XT15.)
2. Describe the current and past problems / symptoms  
\_\_\_\_\_  
\_\_\_\_\_
3. When is the last time you had your instrument serviced? \_\_\_\_\_  
What was done? \_\_\_\_\_
4. What is the fault displayed on the LCD screen? \_\_\_\_\_
5. How often does the fault occur? \_\_\_\_\_
6. How often is the instrument run? \_\_\_\_\_
7. Solvent used: \_\_\_\_\_ What type of samples do you run? \_\_\_\_\_
8. How many years has the instrument been in service? \_\_\_\_\_
9. Do you hydrolyze your samples before extractions are run? \_\_\_\_\_
10. Water supply: Do you use a chiller or laboratory water? \_\_\_\_\_
11. What is the temperature of the water supplying the XT? \_\_\_\_\_ °C
12. What is the air temperature of your lab? \_\_\_\_\_ °C
13. Turn the instrument OFF and the water supply ON. Do you observe water flowing out of the drain line? YES/NO

**PLEASE PERFORM THE FOLLOWING TEST AND SEND THE RESULTS TO SERVICE@ANKOM.COM**

- 1A. (For XT10) Add 350 ml of solvent to the vessel.
- 1B. (For XT15) Make sure sufficient solvent is in the sight glass, if not, add as needed.
2. Turn Instrument on and run a 20 minute extraction.
3. (For XT15) Record the time to fill the vessel \_\_\_\_\_ seconds.
4. Record temperature (LCD screen) and pressure (gauge) readings every minute until "Process Complete" or a fault appears on the LCD screen.

1. _____ °C/ _____ psi	11. _____ °C/ _____ psi	21. _____ °C/ _____ psi
2. _____ °C/ _____ psi	12. _____ °C/ _____ psi	22. _____ °C/ _____ psi
3. _____ °C/ _____ psi	13. _____ °C/ _____ psi	23. _____ °C/ _____ psi
4. _____ °C/ _____ psi	14. _____ °C/ _____ psi	24. _____ °C/ _____ psi
5. _____ °C/ _____ psi	15. _____ °C/ _____ psi	25. _____ °C/ _____ psi
6. _____ °C/ _____ psi	16. _____ °C/ _____ psi	26. _____ °C/ _____ psi
7. _____ °C/ _____ psi	17. _____ °C/ _____ psi	27. _____ °C/ _____ psi
8. _____ °C/ _____ psi	18. _____ °C/ _____ psi	28. _____ °C/ _____ psi
9. _____ °C/ _____ psi	19. _____ °C/ _____ psi	29. _____ °C/ _____ psi
10. _____ °C/ _____ psi	20. _____ °C/ _____ psi	30. _____ °C/ _____ psi

5. When a temperature reaches 90°C, reflux water will start to flow.
  - a. Record flow rate observed in flow gage \_\_\_\_\_
6. During the run is there any solvent dripping into the vent bottle? Circle: YES/NO
7. During the run, what is the temperature of the water coming out of the drain line? \_\_\_\_\_ °C
8. When run is complete, record the amount of solvent in the following,

Vessel \_\_\_\_\_ ml

Teflon Cup \_\_\_\_\_ ml

Vent Bottle \_\_\_\_\_ ml

Site Glass \_\_\_\_\_ ml

*E-MAIL or FAX THE COMPLETED FORM TO SERVICE@ANKOM.COM or 315-986-8091*